

CDIA™ NT-proBNP Immunofluorescence Test Cassette

Cat. No.: DTSJYJ066

Pkg. Size: 25T

Intended use

The CDIA™ NT-proBNP Immunofluorescence Test Cassette is designed for quantitative determination of NT-proBNP in human serum, whole blood or plasma.

General Description

The N-terminal prohormone of brain natriuretic peptide (NT-proBNP or BNPT) is a prohormone with a 76 amino acid N-terminal inactive protein that is cleaved from the molecule to release brain natriuretic peptide. Both BNP and NT-proBNP levels in the blood are used for screening, diagnosis of acute congestive heart failure (CHF) and may be useful to establish prognosis in heart failure, as both markers are typically higher in patients with worse outcome. The plasma concentrations of both BNP and NT-proBNP are also typically increased in patients with asymptomatic or symptomatic left ventricular dysfunction and is associated with coronary artery disease and myocardial ischemia. This kit is intended to use in quantitative detection of NT-proBNP content in human serum, plasma or whole blood.

Principle of the Test

The test uses an anti-human NT-proBNP monoclonal antibody conjugated with fluorescence latex and another anti-human NT-proBNP monoclonal antibody coated on the test line. After the sample has been applied to the test strip, the fluorescence latex-labelled anti-human NT-proBNP monoclonal antibody binds with the NT-proBNP in sample and forms a marked antigen-antibody complex. This complex moves to the test card detection zone by capillary action. Then marked antigen-antibody complex is captured on the test line by the anti-human NT-proBNP monoclonal antibody. The fluorescence intensity of the test line increases in proportion to the amount of NT-proBNP in sample. Then insert test card into the Fluorescent Immunoassay Analyzer FIA7100. According to the ratio of the fluorescence intensity of control region and test region, the concentration of NT-proBNP in sample will be determined and displayed on the screen. .

Reagents and Materials Provides

1. NT-proBNP Test Cassette, 25T
2. Sample diluent buffer

Creative Diagnostics. All rights reserved

45-1 Ramsey Road Shirley, NY 11967, USA

Tel: 631-624-4882 Fax: 1-631-938-8221

E-mail: info@creative-diagnostics.com

www.cd-diatest.com

3. SD card: Calibration curve information card
4. One instruction

Sample preparation

This test can be used for serum, plasma and whole blood samples. The test should be completed within 1 hour after blood collection.

1. The sample should be homogeneous by inversion before testing.
2. Deliver homogeneous serum, plasma (90 µL) or whole blood (130 µL) into the sample diluent buffer and mix well.

Note: Increase the amount of sample if extreme samples make it flow hard.

Assay operation

1. On the main interface of FIA7100, enter testing interface.
2. Read the SD card information.
3. Enter the sample number and other information.
4. Deliver serum, plasma (90 µL) or whole blood (130 µL) into the sample port on the test card.
5. Insert the test card into FIA7100.
6. Reaction time: 15 minutes. The result will be shown on the screen.

Results

Read the result based on the FIA7100.

Limitation

1. As with all diagnostic tests, a definitive clinical diagnosis should not be made based on the result of a single test. The test results should be interpreted considering all other test results and clinical information such as clinical signs and symptoms.
2. It is possible that technical or procedural errors, as well as other interfering substances in the specimen may cause erroneous results.

Specificity

1. Limit of blank: 18 pg/mL.
2. Accuracy: Factory inspection: In the quality control value range.
Type test: The recovery rate should be in the range of 85% to 115%.
3. Precision: Coefficient of Variation CV% ≤10%.
4. Linear range: 50-25000 pg/mL, R ≥ 0.99.

5. Batch difference: The difference of three batches of the kits is not more than 15%.
6. Stability: After the validity period, the kit can also meet the above 1-4 indicators.

Storage

Store the test device at 4 to 30°C. The kit will be valid in 12 months.

Notice for Operations

1. Please do the assay following the instruction, don not touch the membrane of the strip.
2. This strip is used for only once. Please do not use it repeatedly.
3. Blood sample that can be seen by the naked eyes can interfere with the test and lead to erroneous result.
4. Insert the test card into FIA7100 immediately after delivering sample into the sample port.